

Product Data Sheet

Mouse Anti-Biotin, Monoclonal Antibody

Catalog number : AB75001

Product Description

The purified monoclonal antibody, mouse anti-Biotin monoclonal antibody can be used for detecting the biotin-labeled antibodies, proteins and biotin labeled DAN/RNA probes in immunoassay, IHC, in situ hybridization and other detection procedures.

Intended Use

For Research Use or Manufacturing Purpose Only

Immunogen

Biotin

Product Type

Purified Monoclonal Antibody

Host Species

Mouse

Isotype Classification

IgG1

Cross Reactivity

Not known

Protein Concentration

Each lot is determined by A280 nm extinction coefficient factor of 1.36. The antibody concentration may be prepared at 2 - 10 mg/ml in PBS.

Applications

This product is recommended for use in:

ELISA

in situ hybridization

Storage

Store at 2-8 °C or - 20 °C for long period.