

## Product Data Sheet

# Mouse Anti-BosR (BB0647) Monoclonal Antibody

Catalog number : AB76001

---

### Intended Use

For Research Use Only

### Clone Number

Hybridoma cell line ID: 2F9-D2

### Immunogen

BosR (BB0647) protein of *B. burgdorferi*

### Product Type

Purified Monoclonal Antibody

### Host Species

Mouse

### Storage

Store at 2-8 °C or -20 °C for long term

### Protein Concentration

Each lot is determined by A280 nm extinction coefficient factor of 1.36. The antibody concentration may be prepared up to 1 - 10 mg/ml in PBS.

### Applications

This product is recommended for use in:

ELISA

Western Blot