

Product Data Sheet

Mouse Anti-Digoxigenin, Monoclonal Antibody

Catalog number : AB79295

Product Description

The purified monoclonal antibody, mouse anti-Digoxigenin monoclonal antibody can be used for detecting Digoxigenin-labeled DAN/RNA probes in in situ hybridization, ELISA and other immunodetection procedures.

Intended Use

For Research Use or Manufacturing Purpose Only

Immunogen

Digoxigenin

Product Type

Purified Monoclonal Antibody

Host Species

Mouse

Isotype Classification

IgG1

Cross Reactivity

Not known

Protein Concentration

Each lot is determined by A280 nm extinction coefficient factor of 1.36. The antibody concentration may be prepared at 2 - 10 mg/ml in PBS.

Applications

This product is recommended for use in:

ELISA

in situ hybridization

Storage

Store at 2-8 °C or - 20 °C for long period.