

Product Data Sheet

Goat Anti-Human IgG (H+L), Affinity-Purified

Catalog number : GH9200

Product Description

This affinity-purified antibody binds to heavy and light chains on human IgG. It reacts with all human IgG subclasses. Supplied as a concentrated in PBS buffer, this goat anti-human IgG (H+L) antibody does not contain preservatives or any carrier protein, and is available in large volumes for further manufacturing use.

Product No. GH9200-1: 1 mg in PBS.

Product No. GH9200-10: 10 mg in PBS.

Product No. GH9200-C: bulk qty.

Analysis and Test

The quality of this product is assured via ELISA test.

Intended Use

For Research Use and Manufacturing Purpose Only

Storage

Store at -20 °C.