

# Product Data Sheet

## Goat Anti-Human IgA, Affinity-Purified

Catalog number : GH9300

---

### Product Description

The goat anti-human IgA antibody is affinity-isolated from goat antiserum. The purified anti-human IgA antibody binds to the alpha chains on human IgA and has a minimal or none reactivity to the light chains of human immunoglobulins (IgA, IgG, IgM). It reacts with all human IgA subclasses.

Product No. GH9300-1: 1 mg in PBS. .

Product No. GH9300-C: bulk qty.

### Intended Use

For Research Use and Manufacturing Purpose Only

### Storage

Store at -20 °C.

### Analysis and Test

The quality of this product is assured via ELISA and Immunoblotting tests.