

# Product Data Sheet

## Goat Anti-Rabbit IgG (H+L), Affinity-Purified

Catalog number : GR2200

---

### Product Description

This affinity-purified antibody reacts with heavy and light chains on rabbit IgG. Supplied as a concentrated preparation in PBS buffer, this goat anti-rabbit IgG (H+L) antibody does not contain preservatives or any carrier protein, and is available in large volumes for further manufacturing use

### Analysis and Test

The quality of this product is assured via ELISA and Immunoblotting tests.

### Intended Use

For Research Use or Manufacturing Purpose Only

### Storage

Store at -20 °C.